


Optical Lithography

The creation of different surface structures with sub-micron precision and a high resolution direct-write laser lithography method for various applications.



Direct-write laser lithography can be used to fabricate masks for mask-lithography; draw optic waveguides or micro-lenses. It can be used in nanowire research, microfluidics, and metamaterials as well as to make stamps for holograms or any other application that requires high precision, high-resolution microstructures.

Examples made in our facilities include an organic optical modulator, fabrication of micro lenses, and substrates for measurements of nanowire mechanical properties.


The central device of optical lithography workflow is the micro pattern generator 'μPG 101' from 'Heidelberg Instruments'. Working wavelength $\lambda = 375 \text{ nm}$, substrates can be used up to 6" x 6", structures can be made down to 0.6 μm, recording resolution is 200 nm and max device height is 6mm.

Resists used - positive AZ1518, AZ1050, SPR700 and negative SU-8: GM-1075, SU-8: SU8-2002.

Fabrication of micro lenses



All micro lens device preparation processes take place in a laminar flow cabinet. It starts with cleaning the substrate. The substrate is then covered with photoresist using a spin-coater placed in a laminar flow cabinet, then it is gradually heated to evaporate the solvent and then placed in a 'μPG 101' system to illuminate cylindrical structures with a diameter of

about 40 microns. Photoresist is being exposed after two minutes to create the cylindrical structures. After exposition, photoresist is heated to round the cylindrical structures and create micro lenses. After the fabrication process, micro lenses are examined in Zygo NewView 7100 Optical Surface Profiler to verify the quality of micro lenses.



Substrates for mechanical measurements of nanowires

1-dimensional (1D) nanomaterials (nanowires, nanotubes, nanofibers) have a wide range of applications: from nano-electronics to textiles and composite materials. Mechanical properties (strength, durability, etc.) of 1D nanomaterials are important for multiple applications, such as metal nanowires-based flexible electronics. Nano-mechanical properties of 1D nanostructures can be determined, for example, using atomic force microscopy (AFM) by the 3-point bending tests carried out on these 1D nanostructures deposited on a specific substrate surface. Laser direct writing lithography equipment makes customized substrates for mechanical tests of nanowires using a semiconductor silicon wafer as substrate material.



Simplified technological sequence consists of spin-coating of photoresist on a silicon wafer with a thermal oxide; writing (exposure) the designed pattern using photolithography equipment; chemical development of the photoresist and its baking; etching of unprotected silicon oxide (silica mask formation); anisotropic etching of unprotected regions of the silicon surface; photoresist stripping; dissolution of the remaining silicon oxide and wafer rinsing; cutting of silicon wafer in pieces (each piece is a finished substrate).

Optical waveguides

We have established a workflow for the development and testing of active and passive optical waveguide devices. Waveguides are optical elements used for guiding light. We show an electro-optical (EO) modulator which was prepared using two optical lithography steps. The top and side views of the modulator design are as depicted in Fig. A. It comprises an SU-8 waveguide core, electrodes in the plane with the waveguide core and an EO polymer coating. The waveguide device is then tested using custom developed equipment. An optical image of our SU-8 waveguide splitter can be seen in Fig. B. In Fig. C the device operation with excited first mode in the MZI is shown.

